

COMPLEXATION TITRATION CALCULATION

SDSU CHEM 251

SAMPLE PROBLEM

Plot the titration curve for the titration of 25 mL of 8.40 mM $\text{Ni}(\text{NO}_3)_2$ with a 15 mM EDTA titrant. Note that the solution of nickel nitrate is buffered with 12.0 mM benzylamine hydrochloride and 53.6 mM benzylamine.

Determine the pNi^{2+} values after the following volumes of titrant have been added:

- A) 0.00 mL of titrant
- B) 3.75 mL of titrant
- C) 9.35 mL of titrant
- D) 14.00 mL of titrant
- E) 16.50 mL of titrant

BEFORE THE EQUIVALENCE POINT

- Before the equivalence point the metal is in excess in solution.
- Due to the very high formation constants, a back reaction contribution to the total amount of free metal is insignificant.
- Simply determine the concentration of excess metal ion, keeping in mind the dilution that occurs with the titration.

SAMPLE PROBLEM

Plot the titration curve for the titration of 25 mL of 8.40 mM $\text{Ni}(\text{NO}_3)_2$ with a 15 mM EDTA titrant. Note that the solution of nickel nitrate is buffered with 12.0 mM benzylamine hydrochloride and 53.6 mM benzylamine.

Determine the pNi^{2+} values after the following volumes of titrant have been added:

- A) 0.00 mL of titrant**
- B) 3.75 mL of titrant**
- C) 9.35 mL of titrant**
- D) 14.00 mL of titrant
- E) 16.50 mL of titrant

AT THE EQUIVALENCE POINT

- The dissociation of the complex, which varies in proportion to the pH, is the sole contributor to the free metal in solution.
- Determine the total concentration of metal-EDTA complex that has been formed, accounting for the dilution of the solution.
- Use the conditional formation constant to determine the amount of free EDTA and metal cation.

SAMPLE PROBLEM

Plot the titration curve for the titration of 25 mL of 8.40 mM $\text{Ni}(\text{NO}_3)_2$ with a 15 mM EDTA titrant. Note that the solution of nickel nitrate is buffered with 12.0 mM benzylamine hydrochloride and 53.6 mM benzylamine.

Determine the pNi^{2+} values after the following volumes of titrant have been added:

- A) 0.00 mL of titrant
- B) 3.75 mL of titrant
- C) 9.35 mL of titrant
- D) 14.00 mL of titrant**
- E) 16.50 mL of titrant

PAST THE EQUIVALENCE POINT

- Past the equivalence point the excess of EDTA will gradually shift the equilibrium further towards complete complexation of the metal.
- Use the conditional formation constant, along with the diluted concentration of EDTA-metal complex and excess EDTA to determine the amount of free metal.
- The dissociation of the EDTA-metal complex does not need to be considered, due to the very large formation constants.

SAMPLE PROBLEM

Plot the titration curve for the titration of 25 mL of 8.40 mM $\text{Ni}(\text{NO}_3)_2$ with a 15 mM EDTA titrant. Note that the solution of nickel nitrate is buffered with 12.0 mM benzylamine hydrochloride and 53.6 mM benzylamine.

Determine the pNi^{2+} values after the following volumes of titrant have been added:

- A) 0.00 mL of titrant
- B) 3.75 mL of titrant
- C) 9.35 mL of titrant
- D) 14.00 mL of titrant
- E) 16.50 mL of titrant**